

Foundations of Psychophysiology

Part 5: Electromyography

Dr. Laurens R. Krol
Dr. Marius Klug

2024-05-21



NEUROADAPTIVE
HUMAN-COMPUTER
INTERACTION



Brandenburg
University of Technology
Cottbus - Senftenberg

Psychophysiology

Electromyography (EMG)

History

Physiological origins

Measuring EMG

Analysis of EMG

Psychophysiology: Electromyography

History


Psychophysiology: Electromyography

History

The electric organ in fish was discovered to consist of specialized muscle cells by Francesco Redi and Stefano Lorenzini in the 17th century.

Also because no sparks were visible during movement, shock and muscle activity was assumed to be mechanical in nature.

In 1773, sparks from an eel's shock were finally demonstrated by John Walsh.




Psychophysiology: Electromyography

History

It is with particular satisfaction I make to you my first communication, that the effect of the Torpedo appears to be absolutely electrical;

...

A large Torpedo, very liberal of his shocks, being held with both hands by his electric organs above and below, was briskly plunged into water to the depth of a foot, and instantly raised an equal height into air; and was thus continually plunged and raised, as quick as possible, for the space of a minute. In the instant his lower surface touched the water in his descent, he always gave a violent shock, and another still more violent in the instant of quitting the water in his ascent;




Walsh, J. (1773). *Of the Electric Property of the Torpedo. In a Letter from John Walsh, Esq; F. R. S. to Benjamin Franklin, Esq; LL.D., F. R. S., Ac. R. Par. Soc. Ext., &c.* Philosophical Transactions, 63 (1773-1774), pp. 461-480.

Illustration : Risso, A. (1810). *Ichthyologie de Nice*. Paris, FR: F. Schoell.

History

In 1791, Luigi Galvani discovered that frog legs twitched when subjected to electrical current.

Ultimately, Emil du Bois-Reymond (1849) used a galvanometer to demonstrate increased electric activity during muscle contractions.



Psychophysiology: Electromyography

Physiological origins

Muscle types


- Smooth muscle (involuntary, ANS-controlled)
- Cardiac muscle (striated; involuntary, ANS-controlled)
- Skeletal muscle (striated; voluntary, SNS-controlled)

Striated refers its striated appearance under a microscope.




Voluntary control?


- Facial expressions
- Shivering when cold
- Muscle tension due to disorders and stress
- Relaxation during sleep




Skeletal muscle




Sarcomere




Muscle fiber and myofibril




Sarcomere



Sarcomere contraction



Neuromuscular junction



Motor units

A motor unit consists of a single motor neuron, its axon, and all the muscle fibers innervated by it.

Activation of a motor unit follows an all-or-nothing principle.

Activation of a muscle can vary depending on the amount and size of active motor units.

The ratio ranges from 1:3 to 1:3000, allowing fine or gross motor activity.


Motor unit / muscle fiber action potentials

A motor unit action potential (MUAP) is carried by the neuron, triggering a muscle fiber action potential (MFAP) in the postsynaptic muscle fibers.

Muscle fiber action potentials have

- a slightly lower resting potential (-85 mV),
- a slightly longer duration (up to 5 ms), and
- a slower velocity (3 – 5 m/s).

Origin of EMG



Motor unit potentials

A motor unit potential (MUP) refers to the sum of all muscle fiber action potentials of a single motor unit. This is what EMG measures.

The measured signal is thus a function of

- the number of detected active motor units, and
- the rate at which they are firing.


Multiple motor units each firing at around a maximum of 50 Hz, but not all firing at the same time, can add up to much higher frequencies.

This is why muscle activity is often broadband signal.


Psychophysiology: Electromyography

Measuring EMG

The EMG signal: Raw



The EMG signal: Spectrum



Recording

The amplifier needs to be capable of recording broad spectrum signal with most power in 40 – 100 Hz.

Good ground site important to filter out overlapping mains noise.

High-pass filter can be applied to filter out eye activity when recording from the face.

Subdermal electrodes can be used for precise measurement.

Abrade the skin before applying surface electrodes. Smaller sizes can be used for face recordings.

Electrode placement



Usually a bipolar recording with one electrode placed in the middle, the other at the end of the muscle of interest

A few mm can make a big difference.

- Absolute value are therefore mostly meaningless.
- Repeated measures are difficult, even within subjects.
- Focus on within-subject, within-session differences.

Psychophysiology: Electromyography: Measurement

Electrode placement: Face



Left: Cacioppo, J. T. & Tassinari, Louis G. (1990). Principles of psychophysiology: Physical, social, and inferential elements. Cambridge, UK: Cambridge University Press.

Right: Davis, J. F. (1959). Manual of surface electromyography (WADC Tech. Rep. No. 59-184). Dayton, OH, USA: Wright Air Development Center

Facial muscle relevance


Corrugator: frown

Orbicularis: bright light,
averse stimuli

Levator: disgust


Mentalis: sleep

Zygomaticus: smile



Psychophysiology: Electromyography: Measurement


Electrode placement: Body



Left: Cacioppo, J. T. & Tassinari, Louis G. (1990). Principles of psychophysiology: Physical, social, and inferential elements. Cambridge, UK: Cambridge University Press.

Right: Davis, J. F. (1959). Manual of surface electromyography (WADC Tech. Rep. No. 59-184). Dayton, OH, USA: Wright Air Development Center

Body muscle relevance



Splenius (neck):
aversive stimuli

Trapezius: stress

Tibialis (leg): stress

Artefacts


Aside from the obvious, in particular

- posture,
 - movement artefacts,
 - electrode shifts due to muscular activity, and sometimes
 - ECG, depending on electrode positions,
- need to be taken into account.

Psychophysiology: Electromyography


Analysing EMG

The EMG signal: Raw



Psychophysiology: Electromyography: Analysis

EMG processing



1. Raw signal

2. Rectification

- Half-wave: Only retain positive samples
- Full-wave: take absolute values

3. Smoothing

4. Integration





EMG processing: Integration

Integration reflects a measure of “total EMG output over a given period of time”.

The *period of time* is thus a free parameter, dependent on the research question, and can vary from milliseconds to seconds.



Psychophysiology: Electromyography: Analysis

Example




Psychophysiology: Electromyography: Analysis

Example



Psychophysiology: Electromyography: Analysis

Example



moral character: e.g., "Mark slows down to avoid the puddle..." ▲ observed data — model fit

immoral character: e.g., "Mark accelerates through the puddle..." ■ observed data — model fit

Psychophysiology

Electromyography

Simultaneous motor fiber action potentials recorded from surface electrodes result in a broad-spectrum signal.

Through integration over time, a general EMG output signal can be constructed.

The activity of different muscles, in particular those in the face, can be correlated to different psychological states.

Psychophysiology

Part 5: Electromyography



Dr. Laurens R. Krol
krol@b-tu.de